

**Sterilization load record sample for sterilizers that produce a print out (up to four loads a day).**

Date: \_\_\_\_\_

On completion of each cycle, verify parameters, sign print out, and attach to this record.

Load 1	Load 2	Load 3	Load 4
Cycle selected:	Cycle selected:	Cycle selected:	Cycle selected:
Items in load:	Items in load:	Items in load:	Items in load:
Load started by:	Load started by:	Load started by:	Load started by:
Unloaded by:	Unloaded by:	Unloaded by:	Unloaded by:
Internal chemical indicator change: <input type="radio"/> Yes <input type="radio"/> No	Internal chemical indicator change: <input type="radio"/> Yes <input type="radio"/> No	Internal chemical indicator change: <input type="radio"/> Yes <input type="radio"/> No	Internal chemical indicator change: <input type="radio"/> Yes <input type="radio"/> No
External chemical indicator change: <input type="radio"/> Yes <input type="radio"/> No	External chemical indicator change: <input type="radio"/> Yes <input type="radio"/> No	External chemical indicator change: <input type="radio"/> Yes <input type="radio"/> No	External chemical indicator change: <input type="radio"/> Yes <input type="radio"/> No

**Daily biological indicator testing:**

BI test lot #: \_\_\_\_\_ Expiry date: \_\_\_\_\_

BI control lot #: \_\_\_\_\_ Expiry date: \_\_\_\_\_

Growth is read after \_\_\_\_\_ hours

Time and date BI test and control placed in incubator: \_\_\_\_\_ Initial: \_\_\_\_\_

Time and date BI test and control out of incubator: \_\_\_\_\_ Initial: \_\_\_\_\_

Results: BI test:  positive  negative  
 BI control:  positive  negative

**Routine sterilizer cleaning and maintenance is due on:** \_\_\_\_\_